AMEDA PENGUIN[®]

DELUXE NUTRITIONAL WARMERS | CHAUFFE-LAIT DE LUXE
| CALENTADORES NUTRICIONALES DELUXE

EN | Instructions for use
FR | Mode d'emploi
ES | Instrucciones de uso

 Amēda[®]
Penguin[®]
Nutritional Warmer

Single or Four-Well Deluxe Nutritional Warmers
Chauffe-lait de luxe individuel ou quadruple
Calentadores nutricionales de una cesta o cuatro cestas

Welcome to the Ameda Family!

We are excited to help you meet your nutritional warming needs. With the Ameda Deluxe Single-Well or Four-Well Penguin® Nutritional Warmers, you can safely warm feedings and protect vital nutrients.

The Ameda Deluxe Nutritional Warmers work with the patented design Ameda Therma-Liner® which isolates the feeding container from the water while providing an encompassing kangaroo-like warming environment using water to safely transfer heat.

For any questions regarding your Ameda Deluxe Penguin Nutritional Warmer, CustomerCare™ is available during business hours 8am-5pm CST on weekdays at 1.877.99.AMEDA (26332). Have your serial number for warmers and lot number for Therma-Liners when calling. Please check www.ameda.com for up-to-date hours of operation.

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TABLE OF CONTENTS

1. WARNINGS AND CAUTIONS	4-5
2. INTENDED USE	6
3. PRODUCT DESCRIPTION	6
3.1 Button Description	6
4. SET-UP INSTRUCTIONS.....	7
4.1 Preparing the warmer	7
4.2 Customizing room/bed/pod information.....	7
4.3 Modes	7
4.4 Therma-Liner preparation.....	9
5. INSPECTING BEFORE OPERATION.....	10
6. OPERATING INSTRUCTIONS.....	10
6.1 Warm refrigerated or frozen feedings	10
6.2 Thaw frozen milk to cool liquid state	11
7. CLEANING.....	11
8. MAINTENANCE	12
9. SERVICING AND SPARE PARTS.....	12
10. WARRANTY.....	12
11. TECHNICAL DATA, STORAGE AND OPERATIONAL CONDITIONS	13
12. AMEDA ACCESSORIES.....	13

Icon Glossary

	Warning		Temperature Limitation For Transport/Storage		Serial Number
	Caution		Humidity Limitation For Transport/Storage		Catalogue Number
	Important		Temperature Limitation For Operation		Unit of Measure
	Hot Surface		Humidity Limitation For Operation		CUL Listed
	Grounded		Pressure Limitation For Operation/Transport/Storage		Recyclable (Therma-Liners)
	Keep Dry				
	Manufacturer				

1. WARNINGS AND CAUTIONS

IMPORTANT SAFEGUARDS

When using electrical products, the following basic safety precautions should always be followed.

READ ALL INSTRUCTIONS BEFORE USING THIS PRODUCT

DANGER: To reduce the risk of electrocution:

- **DO NOT** place or store where warmer can fall or be pulled into a sink.
- **DO NOT** place or drop warmer into water or other liquid.
- **DO NOT** reach for electrical product that has fallen into water. Immediately unplug from electrical outlet.

The **WARNING** symbol identifies important safety instructions that if not followed can lead to injury or damage to the warmer. *The warning symbols when used with the following words stand for:*

	WARNING Can lead to serious injury or death.
	CAUTION Can lead to minor injury.
	IMPORTANT Can lead to material damage.

WARNING: To avoid fire, electrocution, serious burns or injury:

- **ENSURE** the power cord is easily accessible and removable, in the event of an emergency requiring immediate and total power disconnection.
- **DO NOT** use the equipment in a manner not specified.
- **DO NOT** use AC power source that applies more than 132 V AC between the supply conductors or between supply conductor and ground.
A protective ground connection by way of the grounding conductor in the power cord is required for safe operation.
- **DO NOT** use parts or accessories other than those that are recommended by the manufacturer.
- **USE ONLY** the power cord and connector appropriate for the voltage and plug configuration in your country.
- If the equipment appears damaged or operates abnormally, protection may be impaired. **DO NOT** attempt to continue operation. When in doubt, have the equipment serviced.

WARNING: To avoid fire, electrocution, serious burns or injury:

- **ALWAYS** operate warmer on a flat, horizontal surface.
- Close supervision is necessary if product is used near children or infants.
- Use product for its intended use **ONLY**, as described in this manual.
- Prior to use, always inspect power cord for damage or exposed wires. If damage and/or exposed wires are found, **DO NOT** use and call Ameda CustomerCare at 1.877.99.AMEDA (26332).
- **NEVER** operate electrical warmer
 - if it has a damaged cord or plug
 - if it is not working properly
 - if it has been dropped, damaged, or becomes wet without being properly and fully assembled
- **ALWAYS** keep cord away from heated surface.
- **ALWAYS** ensure the power cord is disconnected before cleaning or servicing.
- **ALWAYS** ensure the voltage of the power cord is compatible to the power source.
- **ALWAYS** ensure the power cord is inserted in a socket-outlet with a protective earth contact for grounding. Intentional AC interruption is prohibited.
- **ALWAYS** ensure servicing is performed by qualified service personnel only using the Manufacturer's Service Manual. There are no user-serviceable parts.

CAUTION: Can lead to minor injury:

- **DO NOT** use Therma-Liner
 - for more than one shift, not to exceed 12 hours
 - on more than one patient, baby, or child

IMPORTANT: Can lead to material damage.

- **NEVER** put electrical product in water or a sterilizer, as you can cause permanent damage to the nutritional warmer.
- **DO NOT** use abrasive cleaners/detergents when cleaning the nutritional warming unit or accessories. CaviWipe™ or SaniWipe™ are recommended for surface cleaning.
- **DO NOT** use sharp objects on the device buttons or wells.
- **DO NOT** use Therma-Liner for more than one shift or 12 hours.
- **DO NOT** overfill or under-fill Therma-Liner.
- There are no user-serviceable parts. Please ensure servicing is performed by qualified service personnel.

SAVE THESE INSTRUCTIONS

2. INTENDED USE

The Ameda Penguin® Nutritional Warmer is intended to warm feedings to a range that includes body temperature.

3. PRODUCT DESCRIPTION Penguin Nutritional Warmer

Your Penguin Nutritional Warmer is customized with settings and warming profiles according to your unique facility needs. In addition, the gentle vibration of the nutritional warmer helps preserve and suspend feeding nutrients and promote even warming.

3.1 BUTTON DESCRIPTION

Navigation Buttons

Top

Scroll Up

Scroll Up

Middle

Confirm

Go Back

Bottom

Scroll Down

Scroll Down


4. SET-UP INSTRUCTIONS

4.1 PREPARING THE WARMER

1. **PLACE** warmer onto a stable, flat, horizontal surface.
2. **PLUG** the power cord into the power source.
3. **CONNECT** the power cord into the back of the warmer device.

4.2 CUSTOMIZING ROOM/BED/POD INFORMATION

1. **FLIP** the power switch into the on position.
2. **PRESS** and **HOLD** the bottom navigation button.


3. **FOLLOW** the on-screen instructions to **ENTER** up to 5 spaces for Room, up to 5 spaces for Bed and up to 5 spaces for Pod:

- +abc/123
- =return to previous entry
- -abc/123
- Enter button=advance to next entry
- On/Off button=edits complete, return to *Standby mode*

4.3 MODES

STANDBY MODE

The Penguin Nutritional Warmer will be in *Standby mode* and cycle through the following screens whenever powered on and not warming:


WARMING MODES

Fridge to Feed

Feed Temp: range of 90-100°F/32-38°C


Average Time: 5-12 minutes depending on volume

Frozen to Feed

Feed Temp: range of 90-100°F/32-38°C


Average Time: 10-18 minutes depending on volume

Thaw



After Thaw Temp: range of 38-45°F/3-7°C

Time: 15-25 minutes depending on volume

READY TO FEED MODE

The feeding is ready when the screen background turns green similar to the following:

The Remaining Hold Time displays a countdown of the remaining hold time.

The warmer emits a ready tone at the beginning of the mode and with periodic

throughout. The warmer will keep the feeding in a safe temperature range for 30

minutes (until turned off, returned to *Standby mode*, or extends to *Overtime mode*).

THAW COMPLETE MODE


When the *Thaw* feature is used, the feeding should be removed immediately to prevent over-warming.

The warmer emits a ready tone.

The heater immediately shuts off, and milk must be refrigerated.

WARNING: Milk must be refrigerated immediately.

The screen background will turn green and the following phrases will appear on the screen:

OVERTIME MODE

After 30 minutes in *Ready to Feed*, the warmer will enter *Overtime*.

Warming and vibration will turn off.

The screen background will turn red and the word "Overtime" will appear in the middle of the screen:

The Overtime screen will remain on until the unit is turned off or is returned to *Standby* by pushing the On/Off button or Enter button.

Note: As the manufacturer, we recommend not using the feeding after the 30 minute hold time, but that feeding should be handled according to the hospital policy.

DISPLAYED ERRORS

If the Penguin Nutritional Warmer encounters an error, the error will be displayed on a red screen with a warning tone. If this should occur, discontinue use of the warmer and call Ameda CustomerCare at 1.877.99.AMEDA (26332).

4.4 THERMA-LINER PREPARATION

1. At the beginning of a shift, take out one Thermo-Liner per user.
2. Using a permanent marker, enter the following information onto the Thermo-Liner:
 - a. Patient Information
 - b. Date/Time
 - c. Information
3. Open the water compartment (blue zipper) and fill to water line.
4. Pinch one end of the blue zipper between your thumb and forefinger and slide down the length of the zipper to close the water compartment.
5. Repeat to ensure a proper seal. This compartment will remain sealed for the duration of the 12 hour shift. **Note: Use care when removing and re-inserting milk container to prevent damage to interior of water compartment.**

5. INSPECTING BEFORE OPERATION

Inspect the product for damage. Do not use if damaged.
Call Ameda CustomerCare at 1.877.99.AMEDA (26332).

6. OPERATING INSTRUCTIONS

1. Flip the power switch into the on position.
2. Press On/Off Button to begin.
3. Scroll to the warming mode of interest:
 - a. *Fridge to Feed*
 - b. *Frozen to Feed*
 - c. *Thaw*
4. Advance by pressing Enter or Confirm (green check)*.

*Note that if a wrong entry was made, you may press the Go Back selection button to go back a screen.

6.1 WARM REFRIGERATED OR FROZEN FEEDINGS

1. Scroll to choose the container type and size of interest.
2. Press Enter.
3. Only when ready to begin warming the feeding, place the feeding container into the inner compartment (orange zipper) of the Therma-Liner and push it to the bottom. Press firmly along the entire orange zipper to encompass the container as much as possible. Use care when removing and re-inserting feeding containers to prevent damage to the interior water compartment.
4. Place the prepared Therma-Liner in the well with the sealed opening at the top.
5. Press Enter*.
6. The device will begin warming, vibrating, and alternating screens between the profile selected and elapsed time.
7. When the feeding has been warmed to a range that includes body temperature, the device will enter *Ready to Feed mode*.
8. Please remove the feeding from the Therma-Liner and set the Therma-Liner aside until next use during the shift.
9. Press On/Off to return the Penguin to *Standby mode*.

*Note that if a wrong entry was made, you may press the Cancel selection button to go back a screen.

6.2 THAW FROZEN MILK TO COOL LIQUID STATE

CAUTION: THE PENGUIN NUTRITIONAL WARMER IS NOT RECOMMENDED FOR THAWING LESS THAN 15 ml of BREAST MILK IN STORAGE BAGS!

1. Scroll to choose the container type and size of interest.
2. Press Enter.
3. Only when ready to begin warming the feeding, place the feeding container into the inner compartment (orange zipper) of the Therma-Liner and push it to the bottom. Press firmly along the entire orange zipper to encompass the container as much as possible. Use care when removing and re-inserting feeding containers to prevent damage to the interior water compartment.
4. Place the prepared Therma-Liner in the well with the sealed opening at the top.
5. Press Enter*.
6. The device will begin warming, vibrating, and alternating screens between the profile selected and elapsed time.
7. When the feeding has been warmed to a liquid cool state, the device will enter *Thaw Complete mode*.
8. Please immediately remove the feeding** from the Therma-Liner and set the Therma-Liner aside until next use during the shift.
9. Press On/Off to return the warmer to *Standby mode*.

*Note that if a wrong entry was made, you may press the Cancel selection button to go back a screen.

**If there are still some frozen particles in the feeding, let it sit until the remainder thaws, or rotate the container gently to thaw the rest of the way. It is better to have a few remaining ice crystals than to overshoot the refrigerated range. Any portion that is meant for immediate feeding may be placed in a bottle or syringe for further warming. Feeding for the future should be placed immediately in the refrigerator for up to 24 hours. After 24 hours, it should be discarded.

7. CLEANING

1. Disconnect Power Cord from the wall receptacle.
2. Prepare reagent wipes (e.g. CaviWipe™ or SaniWipe™)* containing Isopropanol, Ethylene Glycol Monobutyl Ether (2-Butoxyethanol), Diisobutylphenoxyethoxyethyltrimethyl benzylammonium chloride.
3. Using the prepared reagent wipe, clean and remove all surface debris from the device, including the power cord.
4. Repeat as deemed necessary using a newly prepared reagent wipe each time.
5. After the last cleaning wipe is used, allow device to air dry for a minimum of 10 minutes or dry with a soft cloth before returning the device to service.

8. MAINTENANCE

The Penguin Nutritional Warmer requires no preventative maintenance.

9. SERVICING AND SPARE PARTS

All warmer service and repair can be accomplished by contacting Ameda. CustomerCare will issue a Return Material Authorization (RMA) or provide your facility's biomedical department with parts needed to repair the warmer and return it to operation.

10. WARRANTY

THE AMEDA WARRANTY IS NONTRANSFERABLE. THE AMEDA PENGUIN NUTRITIONAL WARMER IS WARRANTED TO THE ORIGINAL OWNER ONLY.

From the date of the purchase of this product, Ameda, Inc. ("Ameda") warrants the Product against defects in material or workmanship for a period of two years. Ameda's sole obligation under this limited warranty shall be to repair or replace, at Ameda's option, any Product that is determined to be defective by Ameda and determined to be covered by this express limited warranty. Repair or replacement under this limited warranty is the sole and exclusive remedy of the First User. PROOF OF PURCHASE IN THE FORM OF A RECEIPTED INVOICE OR BILL OF SALE EVIDENCING THAT THE PRODUCT IS WITHIN THE WARRANTY PERIOD MUST BE PRESENTED TO OBTAIN WARRANTY SERVICE. This limited warranty is extended by Ameda ONLY to the First User and is not assignable or transferable. For warranty service contact Ameda at 1.877.99.AMEDA (26332).

EXCEPT TO THE EXTENT PROHIBITED BY APPLICABLE LAW, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ON THE PRODUCT IS HEREBY DISCLAIMED. AMEDA SHALL NOT BE LIABLE FOR ANY INCIDENTAL, CONSEQUENTIAL, SPECIAL OR PUNITIVE DAMAGES OF ANY KIND FOR BREACH OF THE EXPRESS LIMITED WARRANTY ON THE PRODUCT OR ANY WARRANTY IMPLIED BY OPERATION OF LAW. OTHER THAN THE EXPRESS LIMITED WARRANTY SET FORTH ABOVE, THERE ARE NO OTHER WARRANTIES THAT ACCOMPANY THIS PRODUCT AND ANY ORAL, WRITTEN OR OTHER REPRESENTATION OF ANY KIND IS HEREBY DISCLAIMED.

The user/owner of this equipment shall have sole responsibility and liability for any injury to person or damage to property resulting from: 1) Operation not in accordance with supplied operating instructions, 2) Maintenance not in accordance with authorized maintenance service instructions, 3) Service by anyone other than manufacturer authorized qualified service personnel, 4) Modification of the equipment or accessories and 5) Use of damaged or unauthorized components and accessories.

11. TECHNICAL DATA, STORAGE AND OPERATING CONDITIONS

Penguin Single Well Nutritional Warmer

Power Rating – 120V AC +/- 10%, 60 Hz, 1.25A
Relative Humidity – 80% up to 31° C to 50% to 40° C
Altitude – 2000 meters
Operating Temperature – 5 to 40° C/41 to 104° F
Storage Temperature – -40 to 75°C/-40 to 167°F

Line Fuse – 2 each, type "T" 1.25A, 250 V AC
Mass – 4.09 kg
Size – 22.2 cm H X 20.5 cm W X 25.9 cm D
Pollution degree 2
Installation category H
Indoor use only

Penguin Four Well Nutritional Warmer

Power Rating – 120V AC +/- 10%, 60 Hz, 4A
Relative Humidity – 80% up to 31° C to 50% to 40° C
Altitude – 2000 meters
Operating Temperature – 5 to 40°C/41 to 104° F
Storage Temperature – -40 to 75° C/41 to 104° F
Line Fuse – 2 each, type "T" 5.0A, 250 V AC
Mass – 9.07 kg
Size – 22.2 cm H X 54.1 cm W X 25.9 cm D
Pollution degree 2
Installation category H
Indoor use only

12. AMEDA ACCESSORIES

You may also order additional warmer accessories: Therma-Liner replacements or Perche Shelf through your sales specialist or Ameda CustomerCare at 1.877.99.AMEDA (26332).

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



Ameda, Inc.
485 Half Day Road
Buffalo Grove, IL 60089
1-866-992-6332(AMEDA)

Mothers Choice Products Ltd.
2133-11871 Horseshoe Way
Richmond, BC V7A5H5
1.800.604.6225

Distribution in other countries: For product information or feedback, call your local distributor or location where you purchased the product. For a listing of distributors in your country, please visit www.ameda.com

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The Ameda Logo, Ameda, Penguin, and/or Therma-Liner are registered trademarks of Ameda, Inc. in the United States and other countries.

The Ameda Logo, Ameda, Penguin, and/or Therma-Liner are registered trademarks of Ameda, Inc. in the United States and other countries.

The Ameda Logo, Ameda, Penguin, and/or Therma-Liner are registered trademarks of Ameda, Inc. in the United States and other countries.

Manufactured by:
Ameda, Inc.
485 H If Day Road
Buffalo Grove, IL 60089
USA

©2017 Ameda, Inc.
1200066-0417

 **Amēda**[®]
Penguin[®]
Nutritional Warmer