Ameda Penguin

Deluxe Nutritional Warmers

EN Instructions for Use 1-17

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 **Ameda Penguin®**
Nutritional Warmer

WELCOME TO THE AMEDA FAMILY

We are excited to help you meet your nutritional warming needs. With the Ameda Deluxe Single-Well or Four-Well Penguin® Nutritional Warmers, you can safely warm feedings and protect vital nutrients.

The Ameda Deluxe Nutritional Warmers work with the patented design Ameda Therma-Liner® which isolates the feeding container from the water while providing an encompassing kangaroo-like warming environment using water to safely transfer heat.

For any questions regarding your Ameda Deluxe Penguin® Nutritional Warmer, CustomerCare™ is available during business hours 8am-4:30pm CST on weekdays at 1.866.99.AMEDA (26332). Have your serial number for warmers and lot number for Therma-Liners when calling. Please check www.ameda.com for up-to-date hours of operation.

TABLE OF CONTENTS

1. WARNINGS AND CAUTIONS.....	4-5
2. INTENDED USE.....	6
3. PRODUCT DESCRIPTION.....	6
3.1 Button Description.....	6
4. SET-UP INSTRUCTIONS.....	7
4.1 Preparing the warmer.....	7
4.2 Customizing room/bed/pod information.....	7
4.3 Customizing Container Menu.....	7-8
4.4 Time to Ready Countdown Timer.....	9
4.5 Modes.....	10-11
4.6 Therma-Liner preparation.....	11
5. INSPECTING BEFORE OPERATION.....	12
6. OPERATING INSTRUCTIONS.....	12
6.1 Warm refrigerated or frozen feedings.....	12
6.2 Thaw frozen milk to cool liquid state.....	13
7. CLEANING.....	13
8. MAINTENANCE.....	14
9. SERVICING AND SPARE PARTS.....	14
10. WARRANTY.....	14
11. TECHNICAL DATA, STORAGE AND OPERATING CONDITIONS.....	15
12. AMEDA ACCESSORIES.....	15
13. EMC INFORMATION.....	16-17

Icon Glossary

Warning

Keep Dry

Serial Number

Caution

Manufacturer

Catalogue Number

Important

Temperature Limitation

Unit of Measure

Hot Surface

Humidity Limitation

CUL Listed

Grounded

Pressure Limitation For Operation/
Transport/Storage

Recyclable
(Therma-Liners)

1. WARNINGS AND CAUTIONS

IMPORTANT SAFEGUARDS

When using electrical products, the following basic safety precautions should always be followed.

READ ALL INSTRUCTIONS BEFORE USING THIS PRODUCT

DANGER To reduce the risk of electrocution:

- **DO NOT** place or store where warmer can fall or be pulled into a sink.
- **DO NOT** place or drop warmer into water or other liquid.
- **DO NOT** reach for electrical product that has fallen into water. Immediately unplug from electrical outlet.

The **WARNING** symbol identifies important safety instructions that if not followed can lead to injury or damage to the warmer. The warning symbols when used with the following words stand for:

WARNING Can lead to serious injury or death.

CAUTION Can lead to minor injury.

IMPORTANT Can lead to material damage.

WARNING To avoid fire, electrocution, serious burns or injury.

- **ENSURE** the power cord is easily accessible and removable, in the event of an emergency requiring immediate and total power disconnection.
- **DO NOT** use the equipment in a manner not specified.
- **DO NOT** use AC power source that applies more than 132 VAC between the supply conductors or between supply conductor and ground. A protective ground connection by way of the grounding conductor in the power cord is required for safe operation.
- **DO NOT** use parts or accessories other than those that are recommended by the manufacturer.
- If the equipment appears damaged or operates abnormally, protection may be impaired. **DO NOT** attempt to continue operation. Have the equipment serviced.

WARNING To avoid fire, electrocution, serious burns or injury.

- **ALWAYS** operate warmer on a flat, horizontal surface.
- Close supervision is necessary if product is used near children or infants.
- Use product for its intended use **ONLY**, as described in this manual.
- Prior to use, always inspect power cord for damage or exposed wires. If damage and/or exposed wires are found, **DO NOT** use and call Ameda CustomerCare at 1.866.99.AMEDA (26332).
- **NEVER** operate electrical warmer
 - if it has a damaged cord or plug
 - if it is not working properly
 - if it has been dropped, damaged, or becomes wet without being properly and fully assembled
- **ALWAYS** keep cord away from heated surface.
- **ALWAYS** ensure the power cord is disconnected before cleaning or servicing.
- **ALWAYS** ensure the voltage of the product is compatible to the power source.
- **ALWAYS** ensure the power cord is inserted in a socket-outlet with a protective earth contact for grounding. Intentional AC interruption is prohibited.
- **ALWAYS** ensure servicing is performed by qualified service personnel only using the Manufacturer's Service Manual. There are no user-serviceable parts.

CAUTION Can lead to minor injury.

- **DO NOT** use Therma-Liner
 - for more than one shift, not to exceed 12 hours
 - on more than one patient, baby, or child

IMPORTANT Can lead to material damage.

- **NEVER** put electrical product in water or a sterilizer, as you can cause permanent damage to the nutritional warmer.
- **DO NOT** use abrasive cleaners/detergents when cleaning the nutritional warming unit or accessories. CaviWipe™ or SaniWipe™ are recommended for surface cleaning.
- **DO NOT** use sharp objects on the device buttons or wells.
- **DO NOT** use Therma-Liner for more than one shift or 12 hours.
- **DO NOT** overfill or under-fill Therma-Liner.
- If servicing is required, ensure servicing is performed by qualified service personnel.

SAVE THESE INSTRUCTIONS

2. INTENDED USE

The Ameda Penguin® Nutritional Warmer is intended to warm feedings to a range that includes body temperature. For use in Hospitals or Institutional settings only. Not for home use.

3. PRODUCT DESCRIPTION Penguin Nutritional Warmer

Your Penguin Nutritional Warmer is a feeding warming device and is customized with settings and warming profiles according to your unique facility needs. The nutritional warmer operates with gentle vibration which helps preserve and suspend feeding nutrients and promote even warming.

3.1 BUTTON DESCRIPTION

Navigation Buttons

Top		Scroll Up		Scroll Up	
Middle		Select		Go Back	
Bottom		Scroll Down		Scroll Down	

Beak Buttons


4. SET-UP INSTRUCTIONS

4.1 PREPARING THE WARMER

1. **PLACE** warmer onto a stable, flat, horizontal surface.
2. **PLUG** the power cord into the power source.
3. **CONNECT** the power cord into the back of the warmer device.

4.2 CUSTOMIZING ROOM/BED/POD INFORMATION


1. **FLIP** the rear power switch into the “**ON**” position.
2. **PRESS** and **HOLD** the bottom navigation button.
3. **FOLLOW** the on-screen instructions to **ENTER** up to 5 spaces for Room, up to 5 spaces for Bed and up to 5 spaces for Pod:
 - +abc/123
 - =return to previous entry
 - -abc/123
 - Enter button=advance to next entry
 - On/Off button=edits complete, return to Standby mode

4.3 CUSTOMIZING CONTAINER MENU


The Penguin Nutritional Warmer (software version 3.14 or later) container selection menu can be customized to suit the needs of each user. Container sizes can be added or removed from the menu for ease of navigation and use.

1. **PRESS** and **HOLD** the middle navigation button.
2. This will lead to a **CUSTOMIZATION** screen.
 - a. Scroll to choose the **WARMING MODE** you wish to edit.


- b. Use the **ENTER** button to select the **WARMING MODE** of choice.

3. Scroll to choose the **CONTAINER TYPE** you wish to edit.
 - a. Use the **ENTER** button to select the **CONTAINER TYPE** of choice.


4. Scroll to choose the **CONTAINER SIZE** you wish to enable or disable.
 - a. Press the **ENTER** button to enable or disable the desired container size.
 - Enabled containers are highlighted in **GREEN** (or **BLACK** when selected with the cursor).
 - Disabled containers are shown in the background **PURPLE** color.
 - The middle navigation button can be used to return to a previous section.

Pressing the **ON/OFF** button at any time will return the display to Standby Mode.


To reset all container types and sizes to their enabled state, scroll to **RESET ALL CONTAINERS** on the initial Customization screen and press Enter. A ready tone will sound to confirm the reset.


4.4 TIME TO READY COUNTDOWN TIMER

The Penguin Nutritional Warmer (software version 3.14 or later) features a “Time to Ready” countdown timer for convenience during warming cycles. When enabled, the Time to Ready timer appears on the screen after the system completes the variable heating phase.


1. **PRESS** and **HOLD** the middle navigation button.
2. This will lead to the **CUSTOMIZATION** screen. Scroll to **FEATURES** and press **ENTER**.

3. Press **ENTER** to select **TIME TO READY**

Press **ENTER** until **YES** appears on the screen to enable the timer. Press **ENTER** again to disable the timer, and **NO** will appear.


4. When enabled, the **TIME TO READY** countdown timer will appear on a second line underneath the **ELAPSED TIME** session timer after the system completes the variable heating phase.


4.5 MODES

STANDBY MODE

The Penguin Nutritional Warmer will be in Standby mode and cycle through the following screens whenever powered on and not warming:


WARMING MODES

Fridge to Feed

Feed Temp range of 90-100°F/32-38°C


Average Time 10-20 minutes depending on volume

Frozen to Feed

Feed Temp range of 90-100°F/32-38°C

Average Time 20-45 minutes depending on volume


Thaw

After Thaw Temp range of 38-45°F/3-7°C


Average Time 12-35 minutes depending on volume

READY TO FEED MODE - HOLD TIME (Fridge to Feed and Frozen to Feed only) The feeding is ready when the screen background turns green with the following words displayed:

REMAINING HOLD TIME


Displays a countdown of the remaining hold time. The warmer emits a ready tone at the beginning of the hold time and periodically throughout. The warmer will keep the feeding in a safe temperature range for 30 minutes (until turned off, returned to Standby mode, or extends to Overtime mode).

START TO READY

Displays the total time used for warming the feeding.


THAW COMPLETE MODE

When the Thaw feature is used, the feeding should be removed immediately upon completion to prevent over-warming. The warmer emits a ready tone. The heater immediately shuts off, and milk must be refrigerated. The screen background will turn green and the phrase shown will appear on the screen.

WARNING Milk must be refrigerated immediately.

OVERTIME MODE Fridge to Feed and Frozen to Feed only

After 30 minutes in Ready to Feed HOLD TIME, the warmer will enter Overtime. Warming and vibration will turn off. The screen background will turn red and the word "Overtime" will appear in the middle of the screen.

The Overtime screen will remain on until the unit is turned off or is returned to Standby by PUSHING the On/Off button or Enter button.

Note: As the manufacturer, we recommend not using the feeding after the 30 minute hold time, but that feeding should be handled according to the hospital policy.

DISPLAYED ERRORS

If the Penguin Nutritional Warmer encounters an error, the error will be displayed on a red screen with a warning tone. If this should occur, discontinue use of the warmer and call Ameda CustomerCare at 1.866.99.AMEDA (26332).

4.6 THERMA-LINER PREPARATION

1. At the beginning of a shift, take out one Therma-Liner per user.
2. Using a permanent marker, enter the following information onto the Therma-Liner:
 - a. Patient Information
 - b. Date/Time
 - c. Information
3. Open the water compartment (blue zipper) and fill to water line using cold tap water.
4. Pinch one end of the blue zipper between your thumb and forefinger and slide down the length of the zipper to close the water compartment.
5. Repeat to ensure a proper seal. This compartment will remain sealed for the duration of the 12 hour shift.

Note: Use care when removing and re-inserting milk container to prevent damage to interior of water compartment.

5. INSPECTING BEFORE OPERATION

Inspect the product for damage. Do not use if damaged.
Call Ameda CustomerCare at 1.866.99.AMEDA (26332).

6. OPERATING INSTRUCTIONS

1. Flip the rear power switch into the 'ON' position.
2. Press On/Off button to begin.
3. Scroll to the warming mode of interest:
 - a. Room Temperature/Fridge to Feed
 - b. Frozen to Feed
 - c. Thaw
4. Advance by pressing Enter or Select (green check mark)*.

**Note that if a wrong entry was made, you may press the middle navigation button to go back a screen or press On/Off button to restart.*

6.1 WARM REFRIGERATED OR FROZEN FEEDINGS

1. Scroll to choose the container type and size of interest.
2. Press Enter*.
3. Only when ready to begin warming the feeding, place the feeding container into the inner compartment of the Therma-Liner and gently push it to the bottom. Use care when removing and re-inserting feeding containers to prevent damage to the interior water compartment.
4. Place the prepared Therma-Liner in the well and gently push it to the bottom.
5. Press Enter.
6. The device will begin warming, vibrating, and alternating screens between the profile selected and elapsed time.
7. When the feeding has been warmed to a range that includes body temperature, the device will enter Ready to Feed mode.
8. Please remove the feeding from the Therma-Liner and set the Therma-Liner aside until next use during the shift.
9. Press On/Off to return the Penguin to Standby mode.

**Note that if a wrong entry was made, you may press the middle navigation button to go back a screen or press On/Off button to restart.*

6.2 THAW FROZEN MILK TO COOL LIQUID STATE

Caution: The Penguin Nutritional Warmer is not recommended for thawing less than 15 mL of breast milk in storage bags!

1. Scroll to choose the container type and size of interest.
2. Press Enter*.
3. Only when ready to begin warming the feeding, place the feeding container into the inner compartment of the Therma-Liner and gently push it to the bottom. Use care when removing and re-inserting feeding containers to prevent damage to the interior water compartment.
4. Place the prepared Therma-Liner in the well and gently push it to the bottom.
5. Press Enter.
6. The device will begin warming, vibrating, and alternating screens between the profile selected and elapsed time.
7. When the feeding has been warmed to a liquid cool state, the device will enter Thaw Complete mode.
8. Please immediately remove the feeding** from the Therma-Liner and set the Therma-Liner aside until next use during the shift.
9. Press On/Off to return the warmer to Standby mode.

**Note that if a wrong entry was made, you may press the middle navigation button to go back a screen or press On/Off button to restart.*

***If there are still some frozen particles in the feeding, let it sit until the remainder thaws, or rotate the container gently to thaw the rest of the way. It is better to have a few remaining ice crystals than to overshoot the refrigerated range. Any portion that is meant for immediate feeding may be placed in a bottle or syringe for further warming. Feeding for the future should be placed immediately in the refrigerator for up to 24 hours. After 24 hours, it should be discarded.*

7. CLEANING

1. Disconnect Power Cord from the wall receptacle.
2. Prepare reagent wipes (e.g. CaviWipe™ or SaniWipe™) containing Isopropanol, Ethylene Glycol Monobutyl Ether (2-Butoxyethanol), Diisobutylphenoxyethoxyethylidimethyl benzylammonium chloride.
3. Using the prepared reagent wipe, clean and remove all surface debris from the device, including the power cord.
4. Repeat as deemed necessary using a newly prepared reagent wipe each time.
5. After the last cleaning wipe is used, allow device to air dry for a minimum of 10 minutes or dry with a soft cloth before returning the device to service.

8. MAINTENANCE

The Penguin Nutritional Warmer requires no preventative maintenance.

9. SERVICING AND SPARE PARTS

All warmer service and repair can be accomplished by contacting Ameda. CustomerCare will issue a Return Material Authorization (RMA) or provide your facility's biomedical department with parts needed to repair the warmer and return it to operation.

10. WARRANTY

THE AMEDA WARRANTY IS NONTRANSFERABLE. THE AMEDA PENGUIN NUTRITIONAL WARMER IS WARRANTED TO THE ORIGINAL OWNER ONLY.

From the date of the purchase of this product, Ameda, Inc. ("Ameda") warrants the Product against defects in material or workmanship for a period of two years. Ameda's sole obligation under this limited warranty shall be to repair or replace, at Ameda's option, any Product that is determined to be defective by Ameda and determined to be covered by this express limited warranty. Repair or replacement under this limited warranty is the sole and exclusive remedy of the First User. PROOF OF PURCHASE IN THE FORM OF A RECEIPTED INVOICE OR BILL OF SALE EVIDENCING THAT THE PRODUCT IS WITHIN THE WARRANTY PERIOD MUST BE PRESENTED TO OBTAIN WARRANTY SERVICE. This limited warranty is extended by Ameda ONLY to the First User and is not assignable or transferable. For warranty service contact Ameda at 1.866.99.AMEDA (26332).

EXCEPT TO THE EXTENT PROHIBITED BY APPLICABLE LAW, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ON THE PRODUCT IS HEREBY DISCLAIMED. AMEDA SHALL NOT BE LIABLE FOR ANY INCIDENTAL, CONSEQUENTIAL, SPECIAL OR PUNITIVE DAMAGES OF ANY KIND FOR BREACH OF THE EXPRESS LIMITED WARRANTY ON THE PRODUCT OR ANY WARRANTY IMPLIED BY OPERATION OF LAW. OTHER THAN THE EXPRESS LIMITED WARRANTY SET FORTH ABOVE, THERE ARE NO OTHER WARRANTIES THAT ACCOMPANY THIS PRODUCT AND ANY ORAL, WRITTEN OR OTHER REPRESENTATION OF ANY KIND IS HEREBY DISCLAIMED.

The user/owner of this equipment shall have sole responsibility and liability for any injury to person or damage to property resulting from: 1) Operation not in accordance with supplied operating instructions, 2) Maintenance not in accordance with authorized maintenance service instructions, 3) Service by anyone other than manufacturer authorized qualified service personnel, 4) Modification of the equipment or accessories and 5) Use of damaged or unauthorized components and accessories.

11. TECHNICAL DATA, STORAGE AND OPERATING CONDITIONS

Penguin Single Well Nutritional Warmer

Power Rating: 120VAC +/- 10%, 60 Hz, 1.25A

Relative Humidity: 80% up to 31°C to 50% up to 40°C / 80% up to 87.8°F to 50% up to 104°F

Altitude: 2000 meters

Operating Temperature: 5 to 40°C/41 to 104°F

Storage Temperature: -29 to 70°C/-20 to 158°F

Line Fuse: 2 each, type "T" 1.25A, 250 VAC

Mass: 3.96 kg

Size: 21.8 cm H X 19.4 cm W X 21.9 cm D

Pollution degree 2

Installation category H

Indoor use only

Penguin Four Well Nutritional Warmer

Power Rating: 120VAC +/- 10%, 60 Hz, 4A

Relative Humidity: 80% up to 31°C to 50% up to 40°C / 80% up to 87.8°F to 50% up to 104°F

Altitude: 2000 meters

Operating Temperature: 5 to 40°C/41 to 104°F

Storage Temperature: -29 to 70°C/-20 to 158°F

Line Fuse: 2 each, type "T" 5.0A, 250 VAC

Mass: 9.06 kg

Size: 24.8 cm H X 54.0 cm W X 21.9 cm D

Pollution degree 2

Installation category H

Indoor use only

12. AMEDA ACCESSORIES

You may also order additional warmer accessories: Therma-Liner replacements or Perche Shelf through your sales specialist or Ameda CustomerCare at 1.866.99.AMEDA (26332).

Therma-Liner

Perche Shelf

13. EMC INFORMATION

THIS **ELECTRICAL EQUIPMENT** needs special precautions regarding EMC and needs to be installed and put into service according to the EMC provided in these Instructions for use. Portable and mobile RF communications equipment can affect this **ELECTRICAL EQUIPMENT**. **Warning:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Penguin Nutritional Warmer, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. **Warning:** Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. **Warning:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Guidance and manufacturer's declaration – electromagnetic emissions		
The Ameda Penguin Warmer is intended for use in the electromagnetic environment specified below.		The customer or the user of the Ameda Penguin Warmer should assure that it is used in such an environment.
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Ameda Penguin Warmer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Ameda Penguin Warmer is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity			
The Ameda Penguin Warmer is intended for use in the electromagnetic environment specified below.		The customer or the user of the Ameda Penguin Warmer should assure that it is used in such an environment.	
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for AC line; 100 kHz repetition frequency	± 2 kV for AC line; 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0,5 kV, ± 1kV line-to-line; ± 0,5 kV, ± 1kV, ± 2kV line-to-ground	± 0,5 kV, ± 1kV line-to-line; ± 0,5 kV, ± 1kV, ± 2kV line-to-ground	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<ul style="list-style-type: none"> • 0% Ut; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° • 0% Ut; 1 cycle and 70% Ut; 25/30 cycles Single phase: at 0° • 0% Ut; 250/300 cycles 	<ul style="list-style-type: none"> • 0% Ut; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° • 0% Ut; 1 cycle and 70% Ut; 25 cycles At 0° • 0% Ut; 250 cycles 	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Ameda Penguin Warmer requires continued operation during power mains interruptions, it is recommended that the Ameda Penguin Warmer be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE Ut is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity

The Ameda Penguin Warmer is intended for use in the electromagnetic environment specified below.

The customer or the user of the Ameda Penguin Warmer should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V 0,15 MHz – 80 MHz 80% AM at 1 kHz 6V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80% AM at 1 kHz	3 V 0,15 MHz – 80 MHz 80% AM at 1 kHz 6V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80% AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the Ameda Penguin Warmer, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P}$ 80 MHz to 800 MHz $d = 2,3\sqrt{P}$ 800 MHz to 2,7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz 9 – 28 V/m At IEC 60601-1-2: 2014 Table 9 frequencies	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz 9 – 28 V/m At IEC 60601-1-2: 2014 Table 9 frequencies	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the

electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Ameda Penguin Warmer is used exceeds the applicable RF compliance level above, the Ameda Penguin Warmer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Ameda Penguin Warmer.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Ameda Penguin Warmer

The Ameda Penguin Warmer is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Ameda Penguin Warmer can help prevent electromagnetic interference by maintaining a minimum

distance between portable and mobile RF communications equipment (transmitters) and the Ameda Penguin Warmer as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,7 GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Distributed in Canada by:
Distribue' au Canada par:
Distribuido en Canadá por:

Ameda, Inc.
8206 Cross Park Drive
Austin, TX 78754
1-866-992-6332(AMEDA)

Mothers Choice Products Ltd.
2133-11871 Horseshoe Way
Richmond, BC V7A5H5
1.800.604.6225

Distribution in other countries: For product information or feedback, call your local distributor or location where you purchased the product. For a listing of distributors in your country, please visit www.ameda.com

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The Ameda Logo, Ameda, Ameda Mom Inspired, Hospital Trusted, Flexishield, Finesse™, Cool'N Carry, Store'N Pour, HygieniKit and/or registered trademarks of Ameda, Inc. in the United States and other countries.

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