sports


TENS 2


Instructions for Use
Please read carefully

TensCare™


Quickstart Guide


Connect the device to the electrodes.


Place the electrodes on the area to be treated.


Press and hold the power button to turn the device on.

Select the mode by pressing the MODE button.

Select a programme by pressing the PROG button.

Start the session by increasing the intensity.

Contents

1.	Device Description	4
1.1.	Introduction	4
1.2.	Intended purpose	4
1.3.	Kit contents	4
1.4.	Treatment programmes	4
1.5.	TENS programme	6
1.6.	EMS programme	8
1.7.	Massage programme	13
1.8.	Manual programme	15
1.9.	Controls and display	17
2.	How to Use	18
2.1.	Device Setup	18
2.2.	EMS electrode placement	20
2.3.	Starting a session	25
2.4.	Completing a treatment	26
3.	Troubleshooting	26
4.	Device Care	27
4.1.	Cleaning	27
4.2.	Storage	27
4.3.	Charging Sports TENS 2	27
4.4.	Electrode and lead wire care	29
4.5.	Replacement parts	29
5.	Technical Specification	30

1. Device Description

1.1. Introduction

The **Sports TENS 2** combines the treatment capabilities of a TENS device, an EMS or NMES device, and a relaxing massage stimulator in one unit.

It sends a small and gentle stimulation (similar to your natural nerve impulses) to the muscles being treated via two or four auto-adhesive electrode pads, in order to help you to firm and tone your body.

Rehabilitation: neuromuscular facilitation, muscle re-education, muscle training, prevention/slowing of atrophy/hypotrophy. EMS has been shown to improve blood flow and circulation.

Sports training: improve muscle strength and performance; warm up, strength, speed, power, resistance, endurance and recovery.

Muscle relaxation: loosening up of muscular tension, improving muscular fatigue symptoms and reduced fatigue and muscle tension.

1.2. Intended purpose


Sports TENS 2 is a medical device designed to be used in the home healthcare environment to improve muscle tone and bulk, as well as provide symptomatic relief and management of acute or chronic localised pain and muscle pain. It is suitable for use by all who can control the device and understand the instructions.

Do not use the device for any purpose other than this intended use.

 **Warning:** Not suitable for use in children without medical supervision.

1.3. Kit contents

- 1 x **Sports Tens 2** device
- 4 x 50x50mm electrode pads (E-CM5050)
- 2 x Lead wires (L-ST2)
- 1 x Detachable belt clip
- 1 x Charging adapter
- 1 x Storage pouch
- 1 x Instructions for Use

1.4. Treatment programmes

Each programme has its own combination of Frequency and Pulse Width settings which allow for different muscle groups to be exercised and toned through the electrode pads.

- Frequency (measured in Hz - pulses per second)

For EMS programmes:

Low frequencies (1-10 Hz) have a purifying and relaxing effect through individual contractions.

Medium frequencies (20-50 Hz) can put a high level of strain on the muscle, thus promoting the muscular structure.

Higher frequencies, there is permanent tension in the muscle (tetany) combined with premature muscle fatigue. They are therefore preferably used for elasticity and maximum strength training.

For TENS programmes:

Pain gate: A high frequency of 110 Hz is good at blocking pain signals.

Endorphin release: A low frequency of 4 or 10 Hz allows for the release of endorphins, the body's natural morphine-like substances.

- Pulse Width (measured in μs - millionths of a second)

The **Sports TENS 2** unit has pulse widths of 20 to 350 μs . Generally speaking, the higher the pulse width, the more "aggressive" the stimulation feels, if the pulse width is set high enough, it will usually elicit a muscle contraction, which is required for an effective toning of the muscle group selected.

For the TENS programmes (pain relief), it is not desired to elicit a muscle contraction, however Pulse Rate is important because different frequency settings target different nerve groups and the setting will determine if the "Gate Theory" or "Endorphin Theory" of TENS will be used.

Work/Rest

All the EMS programmes include a WARM UP, TRAIN, and COOL DOWN phase. In the TRAIN phase, the stimulation is intermittent.

Once you have reached a comfortable level, 5 seconds after you stop pressing the button, the intermittent Work/Active Rest phase will start.

For each programme, the unit will go through a stimulation phase (Work) for 1-40 seconds, followed by a phase without stimulation (Active Rest) for a similar time.

This is to allow your muscles to recover from exercise and avoid fatigue.

During the Work period (**W**), the unit will generate contractions, try to contract the muscle voluntarily at the time to get the best possible results.

This cycle Work/Rest will continue for the duration of the programme selected.

W and **R** are displayed on the screen.

Ramp

RAMP is the time in seconds taken to move up and down between zero and the set stimulation strength (i.e. the time taken to move between the Work stimulation and the Rest period). The **Sports TENS 2** has a fixed ramp time of 1.5 s up and 0.75 s down.

Constant and Burst Modes

Constant mode is when the sensation is continuous versus Burst mode when the sensation, as its name implies, is one of on and off. Constant mode is often used for acute pain via Pain Gate effect whereas

Burst mode is useful in chronic pain relief. Burst gives a combination of Pain Gate and Endorphin Release, but the squeezing feeling may not be as comfortable. The stimulation intensity will need to be relatively high. In Constant mode, the sensation is continuous and is more of a tingling pins and needles type.

Modulation Modes

Modulation is when either the frequency or pulse width sweep across range of settings. In modulation mode, the machine delivers a less regular pattern of TENS stimulation in an attempt to reduce or minimise the accommodation effects of regular,

patterned stimulation. This is potentially most useful for patients who use TENS for several hours a day, if for no other reason than accommodation occurs at a slower rate and therefore less intensity adjustment may be required. See programmes F and G.

HAN Stimulation

Invented by Professor Han, it is an alternating mode or shifting between high and low frequencies; stimulation works at a high frequency and a relatively low pulse width; and a low frequency and a relatively high pulse width. This is thought to produce a better analgesic effect and could be used to relieve both acute and chronic pain.

1.5. TENS programme

Prog	Preset / Manual	Frequency (Hz)	Pulse Width (μ s)	Output	Prog time (min)	P or E *
1	Preset	80	150	Constant	C	P
2	Preset	100	200	Constant	C	P
3	Preset	2	250	Constant	30	E
4	Preset	100 (in 2 Hz burst)	150	Burst	30	P+E
5	Preset	150 (in 2 Hz burst)	200	Burst	30	P+E
6	Preset	2/80	200/100	Han	30	P+E
7	Preset	2/100	200/150	Han	30	P+E
8	Preset	10/100	250	PFM	C	P+E
9	Preset	2/120	200/100	FM	C	P+E
10	Preset	100	75	DTENS	C	P

11	Manual	2-150	50-300	Constant	5-90/C	P or E
12	Manua	2-150	50-300	Burst 2Hz	5-90/C	P+E or E

- * P=Pain Gate. Gives relief when being used. Wears off after a few hours.
- E= Endorphin Release. Builds up over 40 mins. Lasts several hours after use.

There is not one programme for a particular condition, and the best choice varies from one person to another, even if they have the same type of pain. Each user needs to select both the programme and the positioning of the electrode pads that is best for them. You may need to try a few positions/programmes before finding the one that suits you.

TENS PRESET PROGRAMMES

Programme 1 & 2:

P – Constant stimulation, the programme generates a uniform smooth sensation. This is the gentlest programme. It is advisable to use for your first trial of TENS while you get used to the sensation.

Programme 3:

E – Pulsing or tapping sensation associated with muscle contractions. Uses the Endorphin release mechanism. It is suitable to be used for chronic pains such as back pain & arthritis, but not for injured muscles where the contractions may delay healing.

This low frequency is also known as Acupuncture-Like TENS and can be used on acupuncture points. One example of this is application on acupuncture point C6 for nausea, travel and morning sickness.

Programme 4 & 5:

P + E – This programme combines P and E, and is thought to be the most effective for radiating pain in arms and legs and deep

muscular pain (Osteoarthritic pain in the knee, Sciatica central pain).

Programme 6 & 7:

P + E – Han stimulation. This programme alternates between low and high frequency every 3 seconds. The effective pain relief is provided by a combination of P and E.

Programme 8:

P+E – PFM, Pulse Width and Frequency Modulation. PW decreases as F increases. This programme gives a pleasant slow relaxing massage sensation. The continuous change helps prevent accommodation occurring.

Programme 9:

P+E – FM, Frequency Modulation, this programme constantly varies the sensation, and is ideal for treating chronic pain over long period (example 5+ hours). Neck, shoulder, Elbow pain, Rheumatic pain, Lumbago, Hip pain, Osteoarthritic pain in the knee.

Programme 10:

P – DTENS: Reduced pulse width is ideal for treating nerve rich areas (tension type headache, facial pain, neck pain, post herpetic neuralgia, migraine).

Programme 11:

Manual programme with constant stimulation.

Programme 12:

Manual programme with burst stimulation

TENS can be used for as long as it is necessary. Continuous treatment is fine, but the electrode pads should be repositioned regularly (at least every 12 hours) to allow the skin to be exposed to the air.

The Endorphin Release programmes work better when the strength is high enough to cause small muscle contractions. Best results are achieved with a session duration between 20 and 40 minutes. Longer use may cause muscle ache.

1.6. EMS programme

PROG	Preset / Manual	Phase	Work						Active Rest				Prog time (min)
			Frequency (Hz)	Pulse width (µs)			Time (s)	Amplitude %	Freq (Hz)	PW (µs)	Time (s)	Amp %	
1	Preset	Warm	6	200	280	340	300	100					6
		Train	10	200	280	340	9	80	4	200	2	50	41
		Cool	3	200	280	240	600	70					10
2	Preset	Warm	6	200	280	340	120	100					2
		Train	22	200	280	340	7	80	4	200	11	70	18
		Cool	3	200	280	340	180	70					3
3	Preset	Warm	6	200	280	340	120	100					2
		Train	25	200	280	340	7	80	4	200	11	50	21
		Cool	3	200	280	340	180	70					3
4	Preset	Warm	6	200	280	340	120	100					2
		Train	40	200	280	340	7	80	4	200	11	50	18
		Cool	3	200	280	340	180	70					3

5	Preset	Warm	6	200	280	340	120	100					2
		Train	55	200	280	340	5	80	4	235	11	50	27
		Cool	3	200	280	340	180	70					3
6	Preset	Warm	5	200	280	340	300	100					5
		Train	60	200	280	340	4	80	5	200	10	50	41
		Cool	3	200	280	340	600	70					10
7	Preset	Warm	5	200	280	340	300	100					5
		Train	65	200	280	340	4	80	4	200	11	25	20
		Cool	3	200	280	340	600	70					10
8	Preset	Warm	6	200	280	340	120	100					
		Train	75	200	280	340	4	80	4	200	14	50	25
		Cool	3	200	280	340	180	70					
9	Preset	Warm	5	200	280	340	300	100					
		Train	100	200	280	340	6	80	4	200	36	25	28
		Cool	3	200	280	340	600	70					
10	Manual	Warm	6	200	280	340	120	100					
		Train	10-120	100-350			1-30	80			1-60*		1-90
		Cool	3	200	280	340	180	70					
10	Default		35	280			5				9		20
11	Manual	Warm	6	200	280	340	120	100					
		Train	10-120	100-350			1-30	80			1-60*		1-90
		Cool	3	200	280	340	180	70					
12	Default		50	300			5						

Before selecting an EMS programme, you will need to select a mode:

MODE I is for small muscles like the face and hands.

MODE II is for medium sized muscles like arms and feet.

Mode III is for large muscles like thighs, buttocks, and abdomen.

EMS PRESET PROGRAMMES

Programme 1: Muscle Calming

Relaxing the muscles as much as possible and to promote the body's natural endorphins to promote pain relief and to improve the blood circulation and provide oxygen into the muscle.

Programme 2: Resume Training 1

To promote the slow twitch fibres to build muscle strength to help reduce muscle atrophy ready for resuming training activities. Used for all type of sports.

Programme 3: Resume Training 2

Progress from programme 2 as tolerance increases.

Programme 4: Resistance 1

Improving and increasing the capacity to develop very high level of muscle force over a long period of time. Improving the efficacy of the oxygen consumption at the muscle level and the capacity to withstand toxin accretion, such as lactic acid. For sports activities that require very high levels of prolonged muscle activity: Rowing, Cycling, Middle distance running.

Programme 5: Resistance 2

Progress from programme 4 as tolerance increases.

Alternative application: Lipolysis

Increasing the flow of blood circulation and modifying the metabolism of the lipocytes. To help stimulate the subcutaneous deposits of fat. To assist reduce or eliminate the Orange Peel effect of the skin surface

Programme 6: Maximum Muscle Contraction

To increase muscle bulk and volume and to improve muscle force. Searching for muscular hypertrophy.

Programme 7: Muscle toning 1

Strengthening the muscles, improving blood circulation and capillary bed density. Ideal for applying to the Thigh, Legs, Bottom and Abdomen.

Programme 8: Muscle toning 2

Similar to programme 7, but adding bulk in preference to endurance

Programme 9: Force Output, Anaerobic activity

Increasing the muscle capacity to a level of instantaneous maximum muscle force, changing muscle force into explosive action. Used for all activities requiring maximum muscle output in a very short space of time, such as Judo, short distance sprinting, throwing the discus or shot-put.

All preset TRAIN phases have Ramp up 1.5 s, Ramp Down 0.75 s.

Most Active Rest phases are at 4 Hz, 200 μ s.

EMS MANUAL PROGRAMMES

Warning: Consult your healthcare professional before altering these settings. Correct settings depend on your muscle tone and exercise goals. Inappropriate settings could cause discomfort, undesired muscle balance, or even muscle injury. Correct settings depend on your muscle tone and exercise goals. Inappropriate settings could cause

discomfort or muscle injury or aggravate Urge incontinence.

Programme 10 & 11:

Manual programmes, you can select the frequency, the pulse width, the work and active rest time.

At first use, the default values are shown.

The Warning triangle is displayed if the Active Rest period is less than:

$\text{Rest Time} = \text{Work Time} * (\text{WorkHz} - 16.66) / 16.66 \text{ seconds.}$

This is because muscle fibres can only activate a limited number of times a minute (about 1000) without becoming fatigued.

EMS PROGRAMME USAGE

Recommended use for all body parts except the face:

Generally, 3 sessions per week, over a period of 4 to 6 weeks, gives good initial results depending on the body shaping, toning, sports training or rehabilitation and recovery needs.

Face (Mode I):

Using the **Sports TENS 2** is a non-invasive option for helping lift and/or tighten the soft tissues of the face, improving facial muscle tone and reducing the signs of ageing.

Recommended use for Face: One treatment session daily, 5 days/week, for 12 weeks.

You may notice some positive effects after just one session, but normally it will take about 6 weeks before you notice a difference.

Take care when adjusting the intensity as facial muscles are particularly sensitive; intensity should be increased gradually, starting with a very low level of stimulation (just above perception) and increasing with care until you reach a good level of stimulation, represented by good muscle activation.

On the face, it is not necessary to reach levels of intensity capable of causing discomfort! Great and significant results are obtained through consistency and patience.

Warning: Use the programme of your choice on mode I and never bring the intensity higher than 40 mA.

Always use small electrode pads (see E-CM25) or the face electrode pads* (see E-CM3050).

*Sold separately.

A minimum distance of 3 cm. must be kept between the electrode pad and the eyeball.

Note: You might notice mild skin redness at the site of the electrode pads after the treatment session. It should fade within 10 to 20 minutes.

Very rarely temporary after-effects, such as mild eyelid flutter, have been reported after facial stimulation. This has always resolved within a few days.

GENERAL ADVICE:

Ideally, use **Sports TENS 2** once every two days for each muscle with a day of rest in between. It is possible to work on pairs of muscles, for example thighs and abdominal muscles, treating one set one day and the other the next day.

Warning: Working on too many muscles at the same time is not recommended.

The object of EMS treatment is to produce powerful muscle contractions.

The strength of the current should be increased to **about three times the level at which you can first feel the tingling**, or to as high as you can stand without causing pain.

Increase the intensity until you feel a very smooth tensing under the pads, there should be little visible movement.

Use a low intensity for the first session to help you to understand how the machine works.

Voluntary muscular activity is more effective than stimulation, and it may improve progress if you combine voluntary contraction with stimulation. For best results, always try to contract voluntarily during the work phases.

During the programme and over the next few days the intensity should be increased gradually to produce marked, but not painful muscle contractions.

After treatment tingling sensations may continue or your skin may feel numb, this is normal.

Once good results have been obtained the sessions may be reduced to once or twice a week until the user feels that optimum body shaping has been achieved.

Note: The stimulation should remain comfortable and relaxing at all times.

1.7. Massage programme

PROG	Prog time (min)	Phase	Phase time (min)	1		2		Output
				Frequency (Hz)	Pulse width(μs)	Freq (Hz)	PW (μs)	
1	30	1	30	85	200	135	100	Burst
2	30	1	2	1	200	15	200	Burst
3	30	1	2	25	200	80	200	Burst
4	44	1	4	1	250	4	250	FM
		2	6	20	250			IM
		3	6	4	250			IM
		4	4	1	250	4	250	FM
		5	6	30	250			IM
		6	6	4	250			IM
		7	4	1	250	8	250	FM
		8	6	40	250			IM
		9	6	4	250			IM
		10	4	1	250	8	250	FM
		11	6	50	250			IM
		12	6	4	250			IM
		13	4	1	250	8	250	FM
5	30	1	30	5	300			Constant
		2	30	8	300			Constant
6	30	1	30	5	200			Constant
		2	30	8	200			Constant
7	30	1	30	5	300			Constant
		2	39	8	300			Constant
8		1	5	250	25	250	250	Burst 1Hz
		2	5					Burst 1.25 Hz
		3	5					Burst 1.42 Hz
		4	5					Burst 1.66 Hz

9	30	1	25	200		2/2
		2	30	200		2/2
		3	40	200		1/1
		4	50	200		1/1
		5	70	200		0.5/1
		6	80	200		0.5/1
10	30	1	1	200		15 pulses
		2	2	200		11 pulses
		3	3	200		18 pulses
		4	4	200		11 pulses
		5	5	200		15 pulses
		6	6	200		20 pulses
		7	9	200		26 pulses
		8	11	200		33 pulses
		9	15	200		30 pulses
		10	25	200		1 pulses
		11	15	200		30 pulses
		12	11	200		33 pulses
		13	9	200		26 pulses
		14	6	200		20 pulses
		15	5	200		15 pulses
		16	4	200		11 pulses
		17	3	200		18 pulses
		18	2	200		11 pulses

IM = Intensity Modulation

FM = Frequency Modulation, this programme constantly varies the sensation.

MESSAGE PRESET PROGRAMMES

Massage programmes can be used to relieve stiffness in the muscles. Simply experiment with the different programmes to find the one that you find the most relaxing.

Programme 1 alternates between two strong squeezing/kneading feelings.

Programme 2 alternates between slow and fast tapping.

Programme 3 alternates between gentle and strong kneading.

Programme 4 has five steps of gradually increasing strength of kneading, alternating with tapping.

Programme 5 and 7 alternates between slow and medium-strong tapping.

Programme 6 is a gentler form of programme 5.

Programme 8 is a pumping sensation at about your own heart rate, that slowly increases in speed over the course of the 20 min treatment. It is especially effective with swollen legs and veins.

Programme 9 is squeezing/kneading in 6 steps of increasing intensity.

Programme 10 generates a tapping sensation gradually speeding up, then slowing down.

1.8. Manual programme

When a programme has manual settings available, **MANUAL** will be displayed. The manual programmes are the TENS programmes **11** and **12** as well as the EMS programmes **10** and **11**.

In the manual programmes, you can select and change the following parameters:

Parameter to Change:	Display Symbol Flashes:
Frequency	Hz
Pulse Width	µs
Treatment Timer (minutes)	min
DATE	0 on the right
HOUR	0 on the left
MEM D	0 on the right
MEM H	0 on the left
And for the EMS programmes only:	
Work Time (seconds)	WORK sec
Rest Time (seconds)	REST sec

1. Select the mode TENS or EMS using the **MODE** button, then select one of the manual programmes pressing **P**.
2. Press the **+MENU** to start changing the parameters.
3. In the TENS programmes, you can change the parameters frequency, pulse width, treatment timer, date, hour, MEM D and MEM H.
In the EMS programmes, you can change the parameters: frequency, pulse width, work time, rest time, treatment timer, date, hour, MEM D and MEM H.

4. To adjust the parameter flashing on the screen, use the + and – buttons.
5. To confirm your selection and go to the next (or previous) parameter press **+MENU** or **-MENU**.
6. Once the changes are complete, press the return button  to confirm. The flashing will stop.

MANUAL PROGRAMME LOCK

Once you are back on the main screen after adjusting the parameters, you can protect the manual settings by pressing and holding  button for 5 seconds.

The key symbol  will start flashing on the screen.

If you try to adjust the parameters, the key symbol  will flash.

Press and hold  for 5 seconds to unlock the settings of the programme.

WARNING TRIANGLE

In the EMS manual programmes a Warning triangle  will be displayed if the Rest period is less than $\text{Rest Time} = \text{Work Time} \times (\text{WorkHz} - 16.66) / 16.66$ seconds


DATE AND TIME SETTING

To set date and time, press **+MENU** and cycle through parameters to **DATE**.

Centre left shows the day of month and


centre right shows the time.

When **DATE** appears and the number on the left is flashing, adjust the day with the + and – buttons and confirm with the return button .

When **HOUR** appears and the number of the right is flashing, adjust the time with the + and – buttons and confirm with the return button .

1.9. Controls and display


2. How to Use

2.1. Device Setup

Connect the lead wires

1. Insert the lead wire plug into the base of the unit.
2. Rotate the body of the plug 90° to lock the lead wire in place. Attaching the electrode pads to the lead wires.

3. Push the pin ends firmly into the pigtail ends of the electrode pads.

Lead wire colour coding.

The ends of the lead wires are coloured red or black. When using the EMS and Massage programmes, one of the electrode pads should be positioned over the motor point.

ON/OFF

1. To turn the unit ON, press the ON/OFF button and hold for 3 to 5 seconds until the display shows.
2. To turn the unit OFF, press the ON/OFF button on the top of the unit and hold for 3 to 5 seconds until the display stops

The unit will turn OFF automatically:

- When the Timer reaches zero,
- If it is left at zero strength for more than 5 minutes.

- The unit will start with the programme you used last.


i Note: Always check unit is OFF before applying or removing electrode pads.

The backlight will turn off 10 seconds after the last button press.

MODE CONTROL

Press the **MODE** button to select a mode.

Display will cycle through TENS, EMS I, EMS II, EMS III and MASSAGE (5 choices). Modes EMS I, II, and III are for use on body areas with small, medium and large muscles.


i Note: The intensity must be reset to zero before **MODE** can be changed, otherwise **PAUSE** is activated.

PROGRAMME CONTROL

The button marked P is the programme control. The **Sports TENS 2** has:

- 10 preset TENS
- 2 manual TENS
- 27 preset EMS (9 preset with 3 modes each)
- 6 manual EMS (2 manual with 3 modes each)
- 10 massage programmes

Each time you press and release the P button, the programme changes and is shown on the LCD. The Hz and μ s settings for each programme are also displayed.


Each time you change the programme, the strength level reverts back to zero. This is a safety feature to alleviate any sudden feeling of a surge, as each programme gives a different sensation.

STRENGTH CONTROLS

The buttons marked Ch1▲, Ch1▼, Ch2▲, Ch2▼ are the strength controls.

The left-hand buttons change strength in the left-hand lead, Ch1, and the right-hand buttons change the right-hand lead, Ch2.

To increase strength in steps of 1 mA, press and release the top part of the channel buttons, Ch1▲ and Ch2▲.

The strength levels are shown on the LCD.

The strength control buttons will not operate until the unit is properly connected to you. If the **Sports TENS 2** detects a disconnection in either channel it will automatically return the strength to zero in that channel and flash '00'.

The unit has 99 levels of strength.

If you hold down the strength controls for 3 to 5 seconds, the strength will start increasing continuously.

You may feel nothing over the first few presses. Continue pressing until the sensation is strong but comfortable. At each phase of the programme the intensity automatically adjusts to give a comfortable stimulation, however further increases during use may be necessary if your body becomes used to the sensation. To decrease the strength, press and release the lower part of the buttons marked Ch1▼ and Ch2▼.

The yellow LED on the output socket indicates that there is an active output which can be high enough to be harmful if not used correctly. The display will remain on for 5 seconds after the plug is removed.

TREATMENT TIMER

The default setting for preset programmes is displayed. You can adjust the treatment timer to any duration between 1 and 90 minutes or set it to continuous.

Low battery

An empty battery symbol will show when you need to charge the battery.

Keypad lock

If you do not press any keys for 10 seconds, the keypad will lock so that only the keys marked Ch1▼ or Ch2▼ will work. This is to avoid accidental changes in setting. To unlock, press either Ch1▼ or Ch2▼ keys.

OPEN CIRCUIT CUTOUT

Note: The electrode pads need to be applied directly to the skin before the intensity increase buttons will operate. If the **Sports TENS 2** is not correctly connected to your body, the strength in the channel which is not properly connected will automatically reset to zero. This is to prevent sudden changes if a broken connection is remade.

PAUSE

Pressing the **MODE** button while a programme is in use stops the stimulation and the timer. The pause symbol ►|| is displayed on the screen. Press the **MODE** button again to resume the programme.

MEMORY

To view the memory, press **+MENU** or **-MENU** to cycle through parameters until **MEM D** is displayed (to view the records on that date),

If no programme has been used for more than 19 minutes, 0 is displayed.

When **MEM D** is displayed, use the **+** and **-** buttons to select the day (shown on the left of the screen).

Having selected the required date, press **+MENU** or **-MENU** to cycle to **MEM H** (to display hour of individual recordings on that date).

Then use the **+** and **-** buttons to view the different recordings (up to 3 per day).

For each recording the hour is shown on the right side, and the values are displayed for **MODE, PROGRAMME**, Time of treatment in minutes (0-99), and the average intensity (0-99) used in each channel.

To clear the memory, press and hold the **+MENU** and **-MENU** buttons together for 5 seconds while in memory mode **MEM D** or **MEM H**.

2.2. EMS electrode placement

The electrode pads must always be used in pairs (two electrode pads on each channel), so that the signal can flow in a circuit.

You may feel more sensation in one pad than the other. This is normal – it depends on where the electrode pads are in relation to your nerves.

Use the square electrode pads (see E-CM5050) for all areas except the face where specific optional electrode pads are required (see ECM25 and E-CM3050).

Warning: Always check that the unit is OFF before attaching or removing electrode pads.

Electrode placement for EMS is important for obtaining the best results.

Place two electrode pads over the bulk of the muscle, with one electrode over the muscle's motor point. The motor point is the area on the skin that is located closest to the motor nerve's entry into the muscle – about 1/3 of the way down the muscle from the spine. Here, it is easiest to trigger a contraction by electrical stimulation.

Experiment by moving the electrode pads across the skin until you locate the point over the muscle that gives the cleanest contraction.


Note: The + electrode pad needs to be plugged in to the red lead wire connector, the – refers to the one for the black lead wire connector.

The easiest way is to apply the electrode pads around/near the source of the pain.


TENS activates the nerves best if it travels along the nerve rather than across it. So, place one pad further from the spine than the source of the pain, and one closer.

If the pain is in, or close to, your spine you can place one pad either side of the spine.


Chest


Triceps


Biceps


Shoulders & Traps


Glutes


Abdominals


Hamstrings

Abdominals


Quadriceps


Face treatment with optional electrode pads

You can use the programme 7 (with MODE I) stimulate the muscles around the face muscles. You will notice the muscles working straight away from the slight twitching. The activation of the muscles stimulates the circulation. This relaxing skin care also contributes to an increase in well-being, making you appear more awake and content. To target specific muscles in the face it is best to use the small circular electrode pads (see E-CM25).


Forehead & eye


Cheeks & eye

Jaw & neck

Face


To stimulate all the muscles on one or either side of the face, it is recommended to use the face electrode pads (see E-CM3050).


FACE – Pad placement with the optional face electrode pads (E-CM3050).

Massage electrode placement


Neck & tension headache

Shoulders & upper back

Circulation & swollen legs

Poor return of blood to the heart is a common problem leading to swollen legs and varicose veins. The programme 8 supports what is known as the venous pump. Waste products are discharged which in turn eases the flow of blood.

Pad advice

- The electrode pads supplied are reusable but for single patient use. The adhesive is a peelable hydrogel (water based).
- In order to obtain the best conductivity through the electrode pads always ensure that they are in good condition and tacky.
- Before use make sure your skin is clean and dry.
- Peel the electrode pads from their protective plastic shield by holding and lifting one corner of the pad and pulling. Do not pull on the pigtail wire of the pad.
- After use always replace the electrode pads on the plastic liner and replace them in the re-sealable plastic bag.
- If the electrode pads dry out, then it is best to buy a replacement pack. In an emergency, it may be possible to restore some of the tackiness of the pad by adding a tiny drop of water on each pad and spreading around. If too much water is added, the electrode pad will become too soft. If that happens then it is suggested in order to try and re-establish some adhesiveness, to place

them sticky side up in a refrigerator for a few hours.

- In very hot weather the gel on the electrode pads may become soft. In such cases place them, still on their plastic liners and in their pouch into a fridge until they return to their normal condition.
- The electrode pads provided are latex-free.
- Replace the electrode pads when they lose their stickiness.
- Poor connection may cause discomfort and skin irritation.
- Storage life of an unopened pack of electrode pads is 2 years. This may be affected by very high temperatures or very low humidity.

2.3. Starting a session

1. Before use make sure your skin is clean and dry.
2. Place the electrode pads on the skin, according to electrode placements.
3. If you are using your **Sports TENS 2** while moving, you can attach the unit to your belt using the belt clip or hang it around your neck using a neck cord (not included). For safety, use thread or wool that will snap easily if the unit snags on something.
4. You can safely lean on the electrode pads. The sensation might change.

Warning: Ensure the **Sports TENS 2** is switched OFF before applying the electrode pads on the skin.

TRAINING SESSION

1. Press and hold the ON/OFF button on the control unit for 3 to 5 seconds to switch the control unit on.
2. When switched on for the first time, the **Sports TENS 2** will automatically select the TENS programme 1. Next time, it will automatically select the programme you were using the last time it was switched off.
3. Select one of the 5 modes (TENS, EMS I, EMS II, EMS III or Massage), and then select a programme. Details in section 8, 9 and 10 will help you identify the best programme to suit you.
4. Identify the correct pads positioning for the selected area and apply the electrode pads as described in section 12.2. for EMS programmes, in section 12.3. for Massage programmes and as described on the electrode pads insert for TENS programmes.
5. Once the electrode pads are correctly positioned, adjust the intensity with the buttons Ch1▼, Ch1▲, Ch2▼ and Ch2▲. The LCD display shows the strength of intensity used for each channel.

Note: If the sensation becomes uncomfortable, reduce the intensity.

2.4. Completing a treatment

When the timer reaches zero, your session is complete, and the unit turns off.

1. Check that the control unit is OFF. If it is not, hold down the OFF button to switch off then remove the electrode pads from your skin by holding the electrode pad itself and gently pulling. Be careful not to pull the wires as this may damage them.
2. Replace the electrode pads to their protective plastic shield and return them to the re-sealable plastic bag.

Note: When removing the electrode pads, DO NOT PULL ON LEAD WIRES.

3. Troubleshooting

If the device is not working correctly, use the following troubleshooting guide.

Problem	Possible causes	Solution
Display will not turn on	Flat batteries	Charge battery.
	Damaged battery	Contact supplier. There is no maintenance or calibration other than replacement of gel pads and battery charging. Do not modify the device as this may reduce safety or effectiveness.
Controls will not work	The keypad is locked.	If LOCK is shown on display, press and hold the ▼ button.
Low battery display	Low batteries	Charge battery.
No sensation	Intensity is not strong enough	Increase the intensity. Most users will feel something at a setting below 20 mA.

Output will not increase above zero	There is a connection error leading to an open circuit. Electrodes are not connected to body. Leads are faulty/damaged.	Ensure machine is attached to your body correctly. Try using the second lead wire. Purchase replacement if necessary.
Sudden change in sensation	If you stop the stimulation and continue a few minutes later, the signal will feel quite a lot stronger.	Always return the intensity to zero after disconnecting the leads or electrode pads.

If you are unable to resolve the issue, please contact our customer service team:

Email: customer.service@tenscare.co.uk

Phone: +44 1372 723434

4. Device Care

4.1. Cleaning

- Clean the device and lead wires at least once a week by wiping with a damp cloth and solution of mild soap and water. Wipe dry.
- Do not immerse the device in water.
- Do not use any other cleaning solution than soap and water.

4.2. Storage

- Storage life of an unopened pack of self-adhesive electrode pads is 4 years. This may be affected by very high temperatures or very low humidity.
- The device has no fixed shelf life.
- Unused batteries have a nominal shelf life of 3 years but will usually last longer than this.

4.3. Charging Sports TENS 2

The **Sports TENS 2** is powered by a built-in rechargeable Li-ion battery. A power adaptor with interchangeable plug is included in the kit. The battery should need charging about once a month. The battery should last at least 15 hours at 50mA 300µs 50Hz. To fully charge the battery it will take about 2 hours.

When the battery is running low, the symbol will show on the screen.

Although the display fades as the batteries run down, the strength of the output does not change until the warning is shown.

Charging the battery

1. Connect the charger to the unit and plug it onto the mains socket.
2. The word BATT and three flashing lines will show on the screen display – that means the battery is being charged.

3. When the battery is charged, the word FULL will display

Note: The device should be charged at least once every 3 months to ensure the battery life is maintained.

Note: The device cannot be used while being charged.

Note: Charging voltage: DC 5V, charging current: 90mA.

Battery Warnings

This product is equipped with a Lithium-ion battery. Failure to follow these instructions could cause the lithium-ion battery to leak acid, become hot, explode or ignite and cause injury and /or damage:

- Do NOT pierce, open, disassemble, or use in a humid and/or corrosive environment.
- Remove batteries from the device if it is unlikely to be used for a long period. Some types of batteries may leak corrosive fluid.
- Do not pierce, open, disassemble, or use in a humid and/or corrosive environment.
- Do not expose to temperatures over 60°C(140F).
- Do not put, store or leave near sources of heat, in direct strong sunlight, in a high temperature location, in a pressurized container or in a microwave oven.
- Do not immerse in water or get wet.
- Do not short-circuit.

- Do not connect the device unless the battery cover is in place.
- If battery leakage occurs and comes in contact with the skin or eyes, wash thoroughly with lots of water and immediately seek medical attention.
- Keep batteries out of the reach of children to prevent them from swallowing them by mistake. If swallowed by child, contact doctor immediately.
- Never attempt to recharge an alkaline battery. Risk of explosion.
- Do not mix old, new or different types of batteries as this may lead to battery leakage or low battery indication.
- Disposal: Always dispose of batteries and device responsibly according to local government guidelines. Do not throw batteries onto a fire. Risk of explosion.

4.4. Electrode and lead wire care

- The electrodes are reusable but for single patient use. The electrodes should provide up to 20 uses. This may vary depending on skin condition, humidity and storage between uses.
- To maximise the life of the electrodes, always ensure the electrodes are replaced on to their plastic liner and stored and sealed in the plastic bag. When the electrodes dry, this may negatively impact the performance of the treatment.

- The adhesive is the electrode is a peelable hydrogel (water based). To restore the tackiness of the electrode, add a tiny drop of water and spread it around. Avoid adding too much, as it will make the electrode too soft. If this occurs, place the electrodes sticky side up in a refrigerator for a few hours to regain adhesiveness.
- Replace the electrodes when they lose their stickiness. Poor connection may cause discomfort and skin irritation.
- To ensure safety and effectiveness we recommend using TensCare electrodes only.
- The life of the lead wires depends greatly on use. Always handle the leads with care. We recommend replacing the lead wires regularly (about every 6 months).

E-CMR60130	Breast electrode pads Pack of 4 (2 pairs)
E-CM50100	Large electrode pads (Pack of 4)
L-ST2	Lead wire
X-MDA534627-1000	Charging adaptor
X-MULTIPA-UK	UK socket
X-MULTIPA-EU	EU socket
X-MULTIPA-USA	US socket

4.5. Replacement parts


The following replacement parts are compatible with the **Sports TENS 2** and may be ordered from TensCare at www.tenscare.co.uk or +44(0) 1372 723434.

E-CM5050	50x50mm electrodes (pack of 4)
E-CM25	Small round electrode pads
E-CM3050	Face electrode pads (Pack of 4 pairs)


5. Technical Specification

Waveform	Symmetrical Rectangular Bi-Phasic
Amplitude (over 1 kOhm load)	90 mA zero to peak +ve in 99 steps +/- 10%
Frequency range	1-120Hz
Pulse width range	5-250µs
Output plug	Fully shielded: touch proof
Channels	Dual channel
Power supply	Built-in lithium-ion battery Model: WEL 703048 -900 mAh/3.7v Manufacturer: We Action Electronics Co., Ltd Ratings: 900 mAh, 3.33Wh, 3.7V
Weight	90g (with belt clip)
Dimensions	114 x 56 x 23 mm
Safety Classification	Internal power source. Designed for continuous use. No special moisture protection.
Service life	>5 years
Applied part classification	Type BF
IP rating	IP22
Environmental Specifications	
Operating conditions	Temperature range: 5 to 40°C Humidity: 15 to 90% RH non-condensing Atmospheric pressure: 700 hPa to 1060 hPa
Storage conditions	Temperature range: -25 to +60°C Humidity: 15 to 75% RH non-condensing Atmospheric pressure: 700 hPa to 1060 hPa

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Thank you for using Sports TENS 2...

At TensCare, we strive to provide the best possible products and services. We value your feedback and continuously work to enhance our offerings. Understanding how our products benefit you is important to us. If you have any insights or suggestions, please get in touch.

Contact us

TensCare Ltd.
9 Blenheim Road,
Epsom, Surrey,
KT19 9BE, U.K.

+44 (0) 1372 723434

www.tenscare.co.uk

customer.service@tenscare.co.uk

Follow us

TensCare Ltd

@tenscareltd

@TensCareLtd

TensCare Limited

EU Declaration of Conformity

Tenscare Ltd herewith declare that the product meets the provisions of REGULATION (EU) 2017/745 which apply to it. The medical device has been assigned to class IIa according to Annex IV of the REGULATION (EU) 2017/745. The product concerned has been designed and manufactured under a quality management system according to Annex IX, of REGULATION (EU) 2017/745. This EU declaration of conformity is issued under the sole responsibility of the manufacturer. No "Common Specification" is applicable.

www.tenscare.co.uk

TensCare Ltd.

9 Blenheim Road, Epsom,
Surrey, KT19 9BE, U.K.

EU DoC

Complies with Medical Device
Regulation (EU) 2017/745.

Advena Ltd.

Tower Business Centre,
2nd Flr, Tower Street, Swatar,
BKR 4013, Malta

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