



AMEDA MINIFIL™ COLOSTRUM CONTAINERS

Ameda MiniFil™ are smaller-sized, ready-to-use colostrum and breast milk collection containers designed for NICU to support the critical early days of breastfeeding. MiniFil™ bottles help eliminate the need to transfer between containers, minimizing any potential loss of “liquid gold” and the risk of exposure to contaminants.

THE MINIFIL™ IS CONVENIENT AND MULTI-PURPOSE:




PUMPING/
COLLECTING



SHORT AND
LONG-TERM
STORAGE



WARMING
AND THAWING



MIXING AND
FORTIFYING



FEEDING



TRANSPORTING



TECHNICAL FEATURES:

- Made from food grade polypropylene
- No phthalates or bisphenol A
- Ready to use, clean from manufacturer
- Fits on any standard breast flange
- Available in convenient packs of 20, or case of 200

AVAILABLE IN TWO SIZES

SKU 800H01 | UPC 810725022474

50mL collection container

SKU 800H02 | UPC 810725023105

80mL collection container