



Ameda CleanEase

MICROWAVABLE STEAM CLEANING BAGS

Ameda CleanEase microwavable steam cleaning bags are a convenient and efficient way to ensure breast pump parts are thoroughly sanitized within minutes. Made with Food Grade plastic, the superior polypropylene CleanEase Steam Bag features three layers to keep the bag safe for up to 20 uses.

Hold Here
Ouvre ici
Tira de aquí

Ensure bag is tightly sealed | Asegúrese de que la bolsa esté bien sellada | S'assurer que le sac est bien scellé

CAUTION: (Precaution) Steam Bag
HOT! CALIENTE! CHAUDE!
Steam Only Clean Bags
Vapeur de vapeur / Pour la Vapeur
Évitez la vapeur / Hot vapour

Provides quick sanitization of breast pump kit, parts and accessories
Proporciona una desinfección rápida del kit de extractor de leche, piezas y accesorios.
Assurez rapidement le kit (le kit), les pièces et les accessoires

BPA Free

20 uses per bag
20 uses per bag
20 utilisations par sac

Leakproof zipper with slider closure for easy hold
Cremallera a prueba de fugas con cierre deslizante para una fácil sujeción
Cremallera a prueba de fugas con cierre deslizante para una fácil sujeción

Expandable bottom helps bag stand upright
El fondo expansible ayuda a que la bolsa se mantenga erguida.
Le fond extensible aide le sac à se redresser

Ameda CleanEase

Microwavable Steam Sanitizing Bag
Bolsas desinfectantes de vapor para microondas
Sacs de désinfection à la vapeur pour micro-ondes

Steam vent for safety, useful to drain excess water after use

Leakproof zipper with slider closure for easy hold

Provides quick sanitization of breast pump kit, parts and accessories

20 uses per bag BPA Free

Extra large sized bag to fit all pump kit parts that need sanitization

Expandable bottom helps bag stand upright

Please use precaution. The major risk with steam bags is scalding, so moms and caregivers should be very careful when they take the bag out of the microwave. || Por favor, tenga precaución. El mayor riesgo con las bolsas de vapor es escaldar, por lo que las madres y los cuidadores deben tener mucho cuidado al sacar la bolsa del microondas. || Si vous plaît utiliser des précautions. Les sacs à vapeur présentent un risque majeur de brûlure; les mères et les soignants doivent donc faire beaucoup attention lorsqu'ils sortent le sac de leur micro-ondes.

CLEANEASE IS CONVENIENT AND MULTI-PURPOSE:

- Sanitize breast pump kit parts such as flanges, flange adapters, diaphragm, valves
- Sanitize Ameda Nipple Shields, DuoShell, collection bottles, locking rings and lids
- Carry breast pump kit parts and accessories while on-the-go

FOR HEALTH AND SAFETY

Ameda meets the standard recommended by HMBANA and the CDC, per the guidelines “Best Practice for Expressing, Storing and Handling Human Milk in Hospitals, Home and Child Care Setting.”

- Major risk with steam bags is scalding, so moms and care givers should be very careful when they take the bag out of the microwave.

AVAILABLE IN 2 SIZES

SKU 800H03 | UPC 810725022030

100 CT bulk pack (for institutions)

SKU 800M03 | UPC 810725022085

7 CT in retail packaging