



## AMEDA MYA JOY DOUBLE ELECTRIC BREAST PUMP


Mya Joy has been designed with efficiency, safety and comfort in mind. It offers the best value, is compact, lightweight and quiet. Mya Joy features the Ameda HygieniKit closed system with the patented diaphragm that protects breastmilk against bacteria, mold and viruses. This pump delivers on the ultimate promise of comfort and ease by offering a wide range of settings, mom-friendly features, and portability.



**PUMP WITHOUT COMPROMISE.**

### KEY FEATURES

- POWERFUL**  
Hospital strength performance with two modes of pumping  
- Stimulation and expression
- PORTABLE**  
Battery-operated, lightweight, can clip right onto a waist band, purse pocket or favorite yoga pants
- SAFE**  
HygieniKit helps protect breast milk from potential contaminants like bacteria, viruses & mold growth
- COMFORTABLE**  
Wide range of suction and speed settings, functions and flange sizes
- DISCREET**  
Operates at less than 45dB and allows for quiet and discreet pumping anywhere





## TECHNICAL FEATURES

- Hospital strength performance with Up to 250 mmHg suction strength
- 2-phase adjustable suction modes: Stimulation (6 Levels) & Expression Mode (12 Levels)
- 2-minute auto switch to Expression Mode
- Large, backlit digital LCD display
- Timer display to easily track pumping time
- Smart auto power off after 45 minutes
- Adaptable for single or double pumping
- Operates on 4 AA batteries or AC adapter
- Features Ameda HygieniKit Milk Collection System
- Compatible with Ameda's CustomFit™ Flange System
- 2-year limited warranty



**131M21**

### AMEDA MYA JOY INCLUDES:

- AC Power Adapter
- Dual HygieniKit® Milk Collection System
- CustomFit™ Flanges (25mm)
- 2 Bottle Holders
- 2 Extra Valves (for HygieniKit)
- 2 Extra Diaphragms (for HygieniKit)
- Quick Start Guide



**131M22**

### AMEDA MYA JOY "PLUS" INCLUDES ALL MYA JOY ITEMS PLUS:

- CustomFit™ Flanges (25mm, AND 28.5mm)
- 4 Extra Valves (for HygieniKit)
- 2 Extra Diaphragms (for HygieniKit)
- Tote Bag
- Accessory Bag (for Parts)
- Manual Breast Pump Handle
- Breast Milk Storage Guidelines