



Ameda Pearl™

HOSPITAL GRADE, MULTI-USER BREAST PUMP

Powerful, Durable Performance in the Hospital and Home

Combines powerful performance inside with a protective shell on the outside. Pearl is designed to build milk supply amongst pump dependent, late pre-term and term moms. Pearl's modern design combines the durability needed for the multi-user/hospital environment with lifestyle features moms expect. Pearl will serve the needs of hospitals and is the ideal multi-user product for milk banks and rental depots. Pearl is perfectly positioned to support continuity of care in the hospital and home.

Ameda Pearl – Modern Design





Ameda Pearl - Proven Hospital Grade Pumping

Clinical study in process.



PROVEN

Over 90% of moms surveyed were highly satisfied with Pearl



COMFORTABLE

Fully customizable expression mode with 72 possible combinations of speed and suction settings



POWERFUL

Hospital Grade / Multi-User Performance with two modes of pumping



DURABLE

Protective housing designed for easy cleaning between users; 3-year limited warranty



PORTABLE

Lightweight with rechargeable battery and integrated handle



WHISPER QUIET

Discreet pumping anywhere



SMART

Vacusense™ Technology ensures consistent vacuum performance in varying altitudes or changing atmospheric conditions



SAFE

Compatible with Ameda's Universal HygieniKit™, the only closed system pump kit proven via study to protect against contaminants like bacteria, viruses and mold growth.

Ameda Pearl - Technical Specifications

Suction Range	Stimulation: 50 - 150mmHg; 6 user-selectable settings Expression: 30 - 250mmHg; 12 user-selectable settings
Cycles per Minute (CPM)	Stimulation: 80 - 120 CPM; 3 preprogrammed settings Expression: 30 - 48 CPM; 6 user-selectable settings
Power Supply	Medical AC Adapter - Output: 15VDC, 2A, 30W Input: 100-240VAC, 50/60Hz, 1A 2600mAh Lithium-ion rechargeable battery; 3 hours per charge
Unit Weight	1.6 Kg
Unit Dimensions	20 x 18 x 15 cm (L x W x H)
Preprogrammed pumping	Mode auto switches from stimulation to expression at 2 mins Auto Power off at 60 mins Vacusense Technology
Warranty	3 year limited



Pearl Breast Pump with AC adaptor

201W01

Pearl Trolley

201H01



For more information
about Ameda Pearl



Ameda Pearl Hard Shell
Transport & Storage Case

501M03



Soft Carry Bag

501M02