



Sterifeed Colostrum Collectors

	Description	Details
	<p>1ml Colostrum Collector</p> <p>Collection of every hand expressed, precious drop of colostrum can be difficult. This 1ml Colostrum Collector offers a tried and tested method of collecting each droplet. Clear markings allow the user to see the volume collected and each collector comes with a stopper so the colostrum can be safely stored. The plunger drives to the end of the barrel using a unique tip to allow the removal of all liquid contained in the syringe, ensuring all of the precious Colostrum can be extracted.</p> <p>Product Specification:</p> <ul style="list-style-type: none"> • Sterile and ready to use • Push on end stopper included • Maximum volume 1ml (graduated) • Manufactured from Polypropylene (BPA Free) • For single use only • Suitable for labelling and freezing 	<p>Order #: 11201</p> <p>Quantity: 100 / Case</p>
	<p>Colostrum Collection and Storage Container</p> <p>Collection of every precious drop of colostrum can be difficult. This collection container has a sufficiently wide opening to allow each droplet to be caught during hand expression. Its specially curved base enables simple and effective removal of colostrum with a colostrum collector or syringe for feeding.</p> <p>Product Specification:</p> <ul style="list-style-type: none"> • Sterile and ready to use • Each sterile package includes 2 containers • Screw on lid • Maximum volume 20ml • Manufactured from Polypropylene (BPA Free) • Box quantity: 50 packs of 2 • Suitable for labelling and freezing • Multi-use container, not just for breastmilk 	<p>Order #: 11200</p> <p>Quantity: 50 / Case</p>