



**SteriFeed**  
Caring for your baby

## STERIFEED IS A MARKET LEADER OF MILK STORAGE BOTTLES

With a large range of high quality products and dedication to research and development, Sterifeed has become a chosen supplier for hospitals and milk banks worldwide. Sterifeed's disposable bottle range represents the ultimate technology in bottles. They are made from food grade, clarified polypropylene, and can be recycled, incinerated or landfilled without any harm to the environment.



### KEY FEATURES

- Ultra-smooth interior for easy removal of all milk droplets
- Airtight & leak resistant for maximum shelf life for breast milk when stored in the refrigerator or freezer
- Standard universal neck fits the collection kit of all breast pumps
- Clear bottles with high quality print makes reading graduations easy
- Clearly marked calibrated and accurate graduations every 1 ml
- Bottles are rigid to prevent force feeding
- Bottles do not crack or shatter when dropped
- Free of Bisphenol-A (BPA) Bisphenol-S (BPS), DEHP, Phthalates
- Suitable for freezing and pasteurising
- Ready to use, clean from manufacturer
- ETO gas is NOT used in sterilization process, per current recommendations
- Available in bulk or wrapped in tamper evident PE bag
- Large size range to meet all hospital requirements, 50ml, 100ml, 130ml, 250ml, 500ml

### AVAILABLE IN THE FOLLOWING SIZES:

50 ml 14100 Sterile Disposable Tamper Evident

100 ml 14103 Sterile Disposable Tamper Evident

130 ml 14101 Sterile Disposable Tamper Evident

250 ml 14102 Sterile Disposable Tamper Evident

130 ml 14001 Bulk Sterile Bottle

250 ml 14002 Bulk Sterile Bottle

500 ml 14093 Reusable Mixing Bottle